
Safety Information and Setup Guide

TANITA™ Muscle Analyzer MC-800

Images in this manual are for illustration purposes only.

Usage Conditions

Temperature Range	: 5 - 35°C
Humidity Range	: 30 - 80% RH (without condensation)
Atmospheric Pressure Range	: 70 - 106 kPa

Storage and Transportation Conditions

Temperature Range	: -10 - 60°C
Humidity Range	: 10 - 90% RH (without condensation)
Atmospheric Pressure Range	: 70 - 106 kPa

To avoid malfunctions, avoid storing the equipment in a place with direct sunlight, significant temperature changes, a risk of dampness, a large amount of dust or a risk of vibration or impact, or in the vicinity of flames.

Please read the guide carefully and keep it available for reference.

EN Table of contents

For your safety	3
For accurate measurement	4
Part names and connection procedure	6
Initial setup	8
Troubleshooting	11
Specifications	11

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For Your Safety

This section explains precautionary measures to be taken to avoid injury to the patients and operators of this equipment, and to prevent damage to property. Please familiarize yourself with this information to ensure safe operation of this equipment. Please note that the MC-800 is not a diagnostic equipment.

Warning

Failure to follow instructions highlighted with this mark could result in injury or damage to property.

Caution

Failure to follow instructions highlighted with this mark could result in death or severe injury.

Prohibited

This mark indicates actions that are prohibited.

Required

This mark indicates instructions that must always be followed.

Warning

Prohibited

No access for people with active implanted cardiac devices.

This equipment must not be used on people with pacemakers or active implantable medical devices.

This equipment passes a weak electrical current through the body which could interfere with and cause the malfunction of electrical medical implants, resulting in serious harm.

Do not handle the plug with wet hands.

This may result in electric shock, fire, or current leakage.

Do not modify or perform internal disassembly of this equipment.

There is a risk of electrocution or injury, and precise analysis cannot be guaranteed.

Do not use multiple adapters.

This may cause fire.

Do not stand on the platform if it is wet.

There is a risk of slipping and falling.

Do not use the equipment if you notice any abnormality before or during use.

If you continue to use the product while it is malfunctioning, injury or accidents may occur.

Do not use the equipment if you have any kind of wound or inflammation on any part of your body that comes into direct contact with the equipment.

Do not touch the power contacts immediately after disconnecting the AC adapter.

There is a risk of electrocution.

Caution

Before assembling, do not turn on or use the equipment. Refer to the "MC-800 Assembly Guide" for detailed assembly instructions.

Prohibited

Do not get the equipment wet.

Avoid using on subjects with allergies to metals.

Allergic reactions may be caused by the chromium used in the electrodes of this equipment.

Do not jump on the equipment.

Do not use the equipment at an angle.

Do not apply excessive force to the equipment.

There is a risk of damage to the equipment and personal injury.

Do not insert fingers into any of the gaps or holes.

Do not place items sensitive to magnetic forces near the equipment.

The magnet of the hand grip holder may corrupt data on devices such as USB memory sticks if these are placed near the equipment.

Assist persons with disabilities.

Another person must assist persons with disabilities who may not be able to take a measurement alone.

Required

This equipment is designated a Class B device (mainly for systems intended to be used in indoor environments) and is CE (EMC) certified, but it may affect devices that are sensitive to electromagnetic waves.

Do not place this equipment within 30 cm of medical electrical equipment.

Do not lean against the equipment.

When unplugging the USB Type-C cable, be sure to hold the connector to avoid breaking the cable.

Use the included AC adapter.

Before measurement, wipe the electrodes and the surrounding area with rubbing alcohol.

Measurements are taken with bare hands and feet, so there is a risk of infection.

Clean the equipment after each use.

Wipe off the equipment if dust accumulates or it becomes dirty.


Stand clear of the measured person during measurement to ensure accuracy.

Continually monitor both the measured person and the equipment for anomalies.

If an anomaly in the measured person or equipment is discovered, take appropriate action, such as stopping the equipment, while ensuring the safety of the measured person.

When moving the equipment, note the following.

1. Unplug the USB Type-C cable and the wireless dongle from the equipment.


2. Remove the tablet from the tablet holder.

3. Hold the equipment as shown in the illustration.

- Move the equipment with the help of two or more people.
- Use a dolly to move the equipment (a dolly is not included).
- To prevent damage to the equipment, hold it by the metal part of the base, and avoid holding the resin parts or the hand grip holder.


! When placing the equipment, lower it slowly and be careful not to pinch your fingers.

Interpretation of analysis results (e.g. evaluation of measurements and formulation of exercise programs based on results) must be performed by a professional.

Weight loss measures and exercise based on self-analysis could be harmful to your health. Always follow the advice of a qualified professional.

To prepare for an emergency, do not place any objects around the power outlet while the equipment is in operation. In the event of an emergency, pull out the power plug from the outlet.

For Accurate Measurement

For people to whom any of the following apply, changes in measurement values should be referred to as a reference.

- People with metal implants.
- People taking medications that induce changes in body water (e.g. diuretics, etc.).
- People who are pregnant, on dialysis or experiencing any swelling.

Environment

Use in a stable location.

There may be errors in the measurement values.

Do not install the equipment in a location subject to strong vibrations.

There may be errors in the measurement values.

For further accuracy and security, position the six adjustable feet on the base until they just make contact with the floor.

Do not use this equipment near other products that emit electromagnetic waves.

There is a risk of interference with some devices such as lighting fixtures, medical equipment, and communication equipment (inverter fluorescent lamps, microwave therapy devices, and mobile phones), causing them to malfunction.

Avoid measuring in multiple locations with greatly differing temperatures.

This may cause inaccurate measurements. Allow the equipment to stand for at least 2 hours before using if it is moved to a new location with a temperature difference of 20°C or more.

Conditions

Avoid measuring after exercise or bathing.

This may cause inaccurate measurements. Take measurements after sufficient rest.

Measure at least 2 hours after waking up.

Measure at least 2 hours after a meal.

Avoid measuring if people are dehydrated or have a fever.

Use the equipment under the same conditions and in the same position as much as possible for accurate tracking of changes.

Readings are affected by the level of hydration and position of the body. Use at the same time of day each time, under the same conditions and in the same body position.

It is recommended to go to the toilet before measurement.

Set the correct amount of clothing.

The amount of clothing you set is automatically subtracted from the measured value to display body weight.

The approximate amount of clothing to set is as follows.

- Long-sleeved top and long pants: 1.0 kg
- Short-sleeved top and shorts: 0.5 kg

*Before taking the measurement, remove any outerwear such as coats and remove any heavy items such as smartphones from your pockets.

Posture

Remove dirt and dust from the palms of your hands and the soles of your feet.

Proper measurements cannot be taken if wearing socks or if the palms or soles of the feet are dirty.

If people have thick or dry stratum corneum on your palms or heels, moisten your palms or heels before taking the measurement.

The fat percentage may be displayed as lower than the actual percentage or errors may occur.


If people's hands and feet are cold, warm them up before taking your measurement.

Correct measurement posture

Stand on the platform so that your feet touch the anti-slip treatment at the front and back of the table.


Improper contact between the feet and electrodes may cause the displayed fat percentage to be lower than the actual percentage, or cause an error.

Hold the hand electrodes so that all fingers are touching.

If the grip is too weak, the measurement will be inaccurate. Be careful when measuring elderly people.

NG: Without all fingers attached to the electrodes.

After gripping the grip, lower your hands.

NG: Bent knees and elbows.

Ensure arms are not touching sides and inner thighs are not touching each other during measurement.

If necessary, place a dry towel between arm and side and/or between thighs.

NG: With skin to skin contact.

Do not move during measurement.

This causes inaccurate measurements.

Periodic Checkpoints

TANITA recommends that each facility conduct periodic checks of the equipment. If you notice any abnormality, do not attempt to repair it yourself but contact us or your dealer.

1. Check the following before use.

- Check the equipment for any damage.
- Set the amount of clothing you are wearing using the app installed on your tablet (page 4).

2. Check the following at least daily.

- Check that the equipment is on a stable and level surface, e.g. on firm flooring, not on a thick carpet.

3. Visually inspect the following at least weekly.

- Inspect the equipment for any contamination.
- Inspect all cables, cords, and connector ends for damage or contamination.
- Inspect all safety-related labeling for legibility.
- Inspect all accessories (electrodes, etc.) for wear or damage.

4. Visually inspect the following at least monthly.

- The screws are not loose.

Update settings, replace items, or call for service as necessary according to the results of the visual inspections. Do not use the equipment if you see any signs of damage.


Equipment that has been damaged must be checked for proper operation by qualified personnel before using again.

Do not wipe the equipment with corrosive chemicals (benzene, cleaner, etc.).

Use a mild detergent to clean the equipment.

Part Names and Connection Procedure

Part names

Supplied Items


- Safety Information and Setup Guide (this guide)
- Assembly Guide
- Attachments for the AC adapter (4)

Type-A

Type-C

Type-O


Type-BF


- AC adapter


- Wireless dongle


- USB Type-C cable [Ⓐ]
(for connecting the AC adapter to the equipment, 3.0 m)

- USB Type-C cable [Ⓑ]
(for supplying power to the tablet, 0.5 m)


- Screws (4)

LED status

Solid white	Initialising after the power is turned on.
Solid blue	After the connection to the tablet is completed.
Blinking blue	While pairing for the first time and after pressing the wireless pairing reset button (page 7).
Slow-blinking blue	While trying to connect to an already paired tablet.
Solid green	Equipment ready (Initialisation and pairing complete).
Blinking pink	When an error occurs (e.g., wireless dongle not connected, measurement error).


Connection of a plug on the top of the column

Do not insert your fingers into the USB Type-C connector as there is a risk of pinching.

Insert the USB Type-C cable B into the USB Type-C connector (page 8). Do not connect anything other than the tablet to this cable.

Connection of plugs the platform

Do not connect anything into any of the connection parts other than the AC adapter or wireless dongle (page 8).

Symbols and their meanings

	EMC Directive (2014/30/EU), RoHS Directive (2011/65/EU)
	EU importer
	For indoor use only
	Input, output
	Manufacturer and date of manufacture
	Rated power input, d.c.
	Rated power output, d.c.
	Serial number
	WEEE- Waste Electrical and Electronic Equipment Directives
	No access for people with active implanted cardiac devices
	Crushing of hands/danger of crushing injuries
	Not to be serviced by users
	No leaning against
	Refer to instruction manual (Safety Information and Setup Guide)
	ON (power)
	OFF (power)


Initial Setup

Adjust the adjustable feet to level the equipment.

For accurate weight measurement, ensure the six adjustable feet are set correctly.




 Do not place the tablet on the tablet holder yet. The tablet could fall and be damaged if the equipment is tilted.

1. Slightly tilt the equipment to one side and close all six adjustable feet.


2. Place the equipment at the installation location and check the following.

- All six adjustable feet are on the ground.

- The equipment is level.


- There is no wobbling.


3. If there are any issues, slightly tilt the equipment and adjust the adjustable feet.

Mount the tablet.



1. Rotate the screw and remove the cover.


2. Connect the USB Type-C cable B (0.5 m) into the output port (page 7) on the top of the column.

3. Bundle the USB Type-C cable behind the tablet holder.

4. Rotate the screw to loosen the upper arm of the tablet holder.


5. Mount the tablet to the tablet holder, adjust and secure the arm, and then tighten the screw.


Supports tablets from 10 to 13 inches.

If this equipment fails due to the tablet, it is not covered under warranty.



! Secure the arm firmly to prevent the tablet from falling.

6. Connect the USB Type-C cable to the tablet.


7. Attach the cover and tighten the screw.

If you want to use a commercially available tablet holder, loosen the screws and remove the included tablet holder. Prepare a tablet holder that is compatible with the VESA standard of 75 mm x 75 mm.


Screw (4)

! Be careful not to drop the tablet holder when removing it.

Plug in the wireless dongle.


Connect the wireless dongle to the USB Type-A connector (page 7).


If the dongle is difficult to insert, turn the two screws and remove the cover.

Turn on the power.


1. Connect the USB Type-C cable A (3.0 m) to the AC adapter.


2. Attach the appropriate attachments.

If you have attached the wrong attachment, press and hold the latch, and then slide the attachment in the opposite direction to remove it.


3. Connect the AC adapter to the platform (page 7) and plug in it to power outlet.

If the cable is difficult to insert, turn the two screws and remove the cover.


4. Turn on the power of the platform.

5. When the power is turned on, the LED on the platform turns solid white (it remains solid white until initialisation is complete).

Turn off the power.

Turn off the power of the platform.

Connect to the tablet.

First time connecting to the tablet

1. After turning on the power and completing initialisation, the LED on the platform will blink blue.
2. Follow the instructions on the screen of the app installed on your tablet to complete pairing.
3. After pairing is complete, the LED on the platform will light up blue for 3 seconds, and then turn green.


If you have previously paired with the tablet

1. After turning on the power and completing initialisation, the LED on the platform will slowly blink blue.
2. After connecting is complete, the LED on the platform will light up blue for 3 seconds, and then turn green.

When connecting to a tablet other than the one previously paired

This equipment cannot be paired with multiple tablets simultaneously. To switch to another tablet, press and hold the wireless pairing reset button for 3 seconds or more.

1. With the equipment turned on, press and hold the wireless pairing reset button for 3 seconds or more.
2. The LED on the platform will blink blue.
3. Follow the instructions on the app screen to complete pairing.
4. After pairing is complete, the LED on the platform will light up blue for 3 seconds, and then turn green.

When purchasing a tablet or updating the OS, check our website in advance for compatible models and OS. Failures or malfunctions caused by external connected devices such as tablets are not covered by the warranty.

Troubleshooting

If the information below does not resolve your issue or if you have any questions, please contact us.

Power doesn't turn on.

Make sure the AC adapter is connected properly.

Unable to pair with tablet.

- Make sure the wireless dongle is connected properly.
- Press and hold the wireless pairing reset button (page 7) for more than 3 seconds to pair it with the tablet again.

The LED blinks pink after powering on.

Check if the wireless dongle is securely connected.

Weight measurements are not stable.

Do not move during measurement.

Impedance measurement error.

Incorrect body composition measurements.

This equipment estimates body composition from measured weight and the resistance to electrical flow (impedance). Weight and impedance fluctuate throughout the day, so measured body composition results may not always be consistent.

However, accurate measurements cannot be made under the following conditions:

- Dry hands and feet
- Holding the hand electrodes weakly
- Skin-to-skin contact (armpits, inner thighs, etc.)

For details, see page 4 and 5.

The included wireless dongle (page 6) complies with Directive 2014/53/EU (RED) and related directives.

The full text of the EU declaration of conformity is available at the following internet address:

<https://tanita.eu/declaration-of-conformity>

Specifications

Model Number		MC-800
Rated Voltage		20 V DC
Impedance Measurement	Measurement System	Multi-frequency 8 electrodes
	Measurement Frequency	5 kHz/50 kHz/250 kHz
	Measurement Current	90 μ A or less
	Electrode Materials	Feet: Glass with ITO coating Handgrips: Nickel-chrome plating
Weight Measurement	Measurement Part	Left side of the body/right side of the body/right arm/left arm/ right leg/left leg/foot to foot
	Measurement Range	75.0–1500.0 Ω (0.1 Ω increments)
	Measurement System	Strain Gauge Load Cell
Interface	Maximum Capacity	300 kg (including preset tare value)
	Scale Interval	0.1 kg
	Measurement Range	2.0–300.0 kg
	Accuracy	\pm 0.2 kg
	Preset Tare Value (clothes weight)	0–10.0 kg (0.1 kg increments)
Usage Conditions	Temperature Range	5–35°C
	Humidity Range	30–80% RH (without condensation)
Product Weight		33 kg
Product Size	Overall	W617× H1195 × D800 mm
	Platform	W418× H55 × D418 mm

*Design and specifications are subject to change without notice.

Disposal

This product is an electronic device. Please dispose of this product appropriately, not as general household waste. Be sure to follow the regulations in your area when disposing of this product.

This product meets the following requirements;
1. EMC Directive (2014/30/EU)
2. RoHS Directive (2011/65/EU)

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TANITA Health Equipment H.K. Ltd.
TANITA Europe B.V.
TANITA (Shanghai) Trading Co., Ltd.
TANITA Europe B.V. Manchester Office
TANITA Corporation of America Inc.